

Human R-Spondin3 protein (Recombinant) (STJP000680)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
Molecular Weight	Recombinant human R-Spondin-3 is a protein consisting of 130 amino acid residue subunits, due to glycosylation migrates as an approximately 24kDa protein on reduced SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	84870
Gene Symbol	RSPQ3
UniProt ID	RSPQ3_HUMAN
Immunogen Region	GLN22-VAL146

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--