

Human HER3 (ERBB3) protein (Recombinant) {Biotin} (STJP000716)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	CHO Cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered PBS pH 7.4.
Molecular Weight	Recombinant Human HER3/ERBB3 is a monomer consisting of 640 amino acid residue, and migrates as an approximately 75 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	2065
Gene Symbol	ERBB3
UniProt ID	ERBB3_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--