

Human Growth hormone receptor (GHR) protein (Recombinant) {Biotin} (STJP000760)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/WB
Host / Source	CHO Cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.2 M filtered 20mM Tris, 50mM NaCl, pH 7.4.
Molecular Weight	Recombinant Human Growth Hormone is a monomer consisting of 240 amino acid residue, due to glycosylation migrates as an approximately 50 kDa protein under reducing conditions in SDS-PAGE. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	2690
Gene Symbol	GHR
UniProt ID	GHR_HUMAN
Immunogen Region	ECD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--