

Human IL-1RN protein (Recombinant) (C-His) (STJP000776)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3557
Gene Symbol	IL1RN
UniProt ID	IL1RA_HUMAN
Immunogen Region	Arg26-Glu177
Immunogen Sequence	RPSGRKSSKMQAFRIWDVNQ KTFYLRNNQLVAGYLQGPV NLEEKIDVVPIEPHALFLGI HGGKMCLSCVKSGDETRLQL EAVNITDLSNRKQDKRFAF IRSDSGPTTSFESACPGWF LCTAMEADQPVSLTNMPDEG VMVTKFYFQEDE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

