

Human MANF protein (Recombinant) (No tag) (STJP000778)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7873
Gene Symbol	MANF
UniProt ID	MANF_HUMAN
Immunogen Region	Leu25-Leu182
Immunogen Sequence	LRPGDCEVCISYLGRFYQDL KDRDVTFSPATIENELIKFC REARGKENRLCYYGATDDA ATKIINEVSKPLAHHIPVEK ICELKKKDSQICELKYDKQ IDLSTVDLKKLRVKELKKIL DDWGETCKGCAEKSDYIRKI NELMPKYAPKAASARTDL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

