

Human IL-18 protein (Recombinant) (No tag) (STJP000779)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3606
Gene Symbol	IL18
UniProt ID	IL18_HUMAN
Immunogen Region	Tyr37-Asp193
Immunogen Sequence	YFGKLESKLSVIRNLNDQVL FIDQGNRPLFEDMTDSDCRD NAPRTIFIISMYKDSQPRGM AVTISVKCEKISTLSCENKI ISFKEMNPPDNIKDTKSDII FFQRSVPGHDNKMQFESSY EGYFLACEKERDLFKLILKK EDELGDRSIMFTVQNE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

