

Human HMGB1 protein (Active Recombinant) (C-His) (STJP000784)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M unfiltered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3146
Gene Symbol	HMGB1
UniProt ID	HMGB1_HUMAN
Immunogen Region	Met1-Glu215
Immunogen Sequence	MGKGDPPKPRGKMSSYAFFV QTCREEHKKHPDASVNFSE FSKKCSERWKTMSAKEKGKF EDMAKADKARYEREMKTYIP PKGETKKKFKDPNAPKRPPS AFFLFCSEYRPKIKGEHPGL SIGDVAKKLGEMWNNTAADD KQPYEKKAALKKEYEKDIA AYRAKGKPDAAKKGVVKA EK SKKKKEEEDEEDEDEDEDE EDEEDEDEEDDDDE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

