

## Human UBE2D3 protein (Recombinant) (C-His) (STJP000786)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Formulation         | Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.                                   |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |   |
|--------------------|---|
| Gene ID            | <a href="#">7323</a>  |
| Gene Symbol        | <a href="#">UBE2D3</a>  |
| UniProt ID         | <a href="#">UB2D3_HUMAN</a>   |
| Immunogen Region   | Met1-Met147   |
| Immunogen Sequence | MALKRINKELSDLARDPPAQ CSAGPVGDDMFHWQATIMGP NDSPYQGGVFFLTIHFPTDY<br>PFKPPKVAFITRIYHPNINS NGSICLDILRSQWSPALTIS KVLLSICLLCDPNPDDPLV<br>PEIARIYKTDKYNRISRE WTQKYAM |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|

