

Human VEGF-A/VEGF121 (K321N) protein (Recombinant) (No tag) (STJP000791)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7422
Gene Symbol	VEGFA
UniProt ID	VEGFA_HUMAN
Immunogen Region	Ala207-Arg327
Immunogen Sequence	APMAEGGGQNHHEVVKFMDV YQRSYCHPIETLVDIFQEYP DEIEYIFKPSCVPLMRCGGC CNDEGLECVPTESNITMQI MRIKPHQGQHIGEMSFLQHN KCECRPKKDRARQENCDKPR R

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

