

## Human IL-33 protein (Recombinant) (No tag) (STJP000793)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">90865</a>
Gene Symbol	<a href="#">IL33</a>
UniProt ID	<a href="#">IL33_HUMAN</a>
Immunogen Region	Ser112-Thr270
Immunogen Sequence	SITGISPIEYLASLSTYND QSITFALEDESYEIYVEDLK KDEKKDKVLLSYYESQHPNS ESGDGVDGKMLMVTLSPTKD FWLHANNKEHSVELHKCEKP LPDQAFFVLHNMHSNCVSFE CKTDPGVFIGVKDNHLALIK VDSSENLCTENILFKLSET

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

