

Human IL-36 beta/IL-1F8 protein (Recombinant) (C-His) (STJP000794)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27177
Gene Symbol	IL36B
UniProt ID	IL36B_HUMAN
Immunogen Region	Arg5-Glu157
Immunogen Sequence	REAAPKSYAIRDSRQMVVWL SGNSLIAAPLSRSIKPVT LH LIACRDTEFSDKEKGNMVYL GIKGDLDLCLFCAEIQQKPTL QLKEKNIMDL YVEKKAQKPF LFFHNKEGSTSVFQSVSYPG WFIATSTTSGQPIFLTKERG ITNNTNFYLD SVE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

