

Human IDO-1 protein (Recombinant) (No tag) (STJP000797)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 50mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3620
Gene Symbol	IDO1
UniProt ID	I23O1_HUMAN
Immunogen Region	Ala2-Gly403
Immunogen Sequence	AHAMENSWTISKEYHIDEEV GFALPNPQENLPDFYNDWMF IAKHLPDLIESGQLRERVEK LNMLSIDHLTDHKSQRLARL VLGCIITMAYVWGKGHGDVRK VLPRNIAVPYCQLSKKLELP PILVYADCVLANWKKKDPNK PLTYENMDVLFSDGDCSK GFFLVSLVVEIAAASAIVI PTVFKAMQMQRDITLLKALL EIASCLEKALQVFHQIHDHV NPKAFFSVLRIYLSGWKGN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

