

Human UBE2L3 protein (Recombinant) (C-His) (STJP000799)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied in 50mM HEPES, 200mM NaCl, 10%Glycerol, 1mM TCEP, pH 7.0
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7332
Gene Symbol	UBE2L3
UniProt ID	UB2L3_HUMAN
Immunogen Region	Met1-Asp154
Immunogen Sequence	MAASRRLMKELEEIRKCGMK NFRNIQVDEANLLTWQGLIV PDNPPYDKGAFRIEINFPAE YPFKPPKITFKTKIYHPNID EKGQVCLPVISAENWKPATK TDQVIQSLIALVNDPQPEHP LRADLAEYKDRKKKCKNA EEFTKKYGEKRPVD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------

