

Human IL-36 alpha/IL-1F6 protein (Recombinant) (C-His) (STJP000802)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27179
Gene Symbol	IL36A
UniProt ID	IL36A_HUMAN
Immunogen Region	Met1-Phe158
Immunogen Sequence	MEKALKIDTPQQGSIQDINH RVWVLQDQTLIAVPRKDRMS PVTIALISCRHVETLEKDRG NPIYLGNGNLNLCIMCAKVG DQPTLQLKEKDIMDLYNQPE PVKSFLFYHSQSGRNSTFES VAFPGWFIASSEGGCPLIL TQELGKANTTDFGLTMLF

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

