

Human IL-36Ra/IL-1F5 protein (Recombinant) (C-His) (STJP000805)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26525
Gene Symbol	IL36RN
UniProt ID	I36RA_HUMAN
Immunogen Region	Val2-Asp155
Immunogen Sequence	VLSGALCFRMKDSALKVLYL HNNQLLAGGLHAGKVIKGEE ISVVPNRWLDASLSPVILGV QGGSQLSCGVGQEPTLTLE PVNIMELYLGAKESKSFTFY RRDMLTSSSFESAAYPGWFL CTVPEADQPVRLTQLPENGG WNAPITDFYFQQCD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

