

Human Carbonic anhydrase 2 protein (Active Recombinant) (No tag) (STJP000809)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	760
Gene Symbol	CA2
UniProt ID	CAH2_HUMAN
Immunogen Region	Ser2-Lys260
Immunogen Sequence	SHHWGYGKHNGPEHWHKDFP IAKGERQSPVDIDTHTAKYD PSLKPLSVSYDQATSLRILN NGHAFNVEFDDSDQKAVLKG GPLDGTYRLIQFHFHWGSLD GQGSEHTVDKKKYAAELHLV HWNTKYGDFGKAVQQPDGLA VLGIFLKVGSAPGLQKVVD VLDSIKTKGKSADFTNFDPR GLLPESLDYWTYPGSLTTPP LLECVTWIVLKEPISVSSEQ VLKFRKLNFNNGEGEPEELM

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

