

Human IL36 gamma/IL-1F9 protein (Recombinant) (No tag) (STJP000810)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56300
Gene Symbol	IL36G
UniProt ID	IL36G_HUMAN
Immunogen Region	Ser18-Asp169
Immunogen Sequence	SMCKPITGTINDLNQQVWTL QGQNLVAVPRSDSVTPVTVA VITCKYPEALEQGRGDP IYL GIQNPEMCLYCEKVGEQPTL QLKEQKIMDLYGQPEPVKPF LFYRAKTGRTSTLESVAFPD WFIASSKRDQPIILTSELGK SYNTAFELNIND

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

