

Human Thioredoxin/SASP/TXN protein (Recombinant) (C-His) (STJP000811)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7295
Gene Symbol	TXN
UniProt ID	THIO_HUMAN
Immunogen Region	Val2-Val105
Immunogen Sequence	VKQIESKTAFQEALDAAGDK LVVVDFSATWCGPCKMIKPF FHSLSSEKYSNVIFLEVDVDD CQDVASECEVKCMPTFQFFK KGQKVGFEFSGANKEKLEATI NELV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

