

## Human IL-16 protein (Recombinant) (C-His) (STJP000813)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3603</a>
Gene Symbol	<a href="#">IL16</a>
UniProt ID	<a href="#">IL16_HUMAN</a>
Immunogen Region	Pro2-Ser130
Immunogen Sequence	PDLNSSTDSAASASAASDVS VESTAEATVCTVTLEKMSAG LGFSLEGGKGSLLHGDKPLTI NRIFKGAASEQSETVQPGDE ILQLGGTAMQGLTRFEAWNI IKALPDGPTVIVIRKSLQS KETTAAGDS

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

