

Human CNTF protein (Active Recombinant) (C-His) (STJP000814)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1270
Gene Symbol	CNTF
UniProt ID	CNTF_HUMAN
Immunogen Region	Met1-Met200
Immunogen Sequence	MAFTEHSPLTPHRRDLCSRS IWLARKIRSDLTALTESYVK HQGLNKNINLDSADGMPVAS TDQWSELTEAERLQENLQAY RTFHVLLARLLEDQQVHFTP TEGDFHQAIHTLLLQVAAFA YQIEELMILLEYKIPRNEAD GMPINVGDGGLFEKKLWGLK VLQELSQWTVRSIHDLRFIS SHQTGIPARGSHYIANNKKM

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

