

Human cAspAT/AATC/GOT1 protein (Recombinant) (No tag) (STJP000815)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2805
Gene Symbol	GOT1
UniProt ID	AATC_HUMAN
Immunogen Region	Ala2-Gln413
Immunogen Sequence	APPSVFAEVPQAQPVLVFKL TADFREDPDPRKVNLGVGAY RTDDCHPWVLPVVKKVEQKI ANDNSLNHEYLPILGLAEFR SCASRLALGDDSPALKEKRV GGVQSLGGTGALRIGADFLA RWYNGTNNKNTVPVYSSPTW ENHNAVFSAAGFKDIRSYRY WDAEKRGDLQGFLNDLENA PEFSIVVLHACAHNPTGIDP TPEQWKQIASVMKHRFLFPF FDSAYQGFASGNLERDAWA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

