

Human HTRA2 protein (Recombinant) (No tag) (STJP000816)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27429
Gene Symbol	HTRA2
UniProt ID	HTRA2_HUMAN
Immunogen Region	Ala134-Glu458
Immunogen Sequence	AVPSPPPASPRSQYNFIAADV VEKTAPAVVYIEILDRHPFL GREVPISNGSGFVVAADGLI VTNAHVVADRRRVVRLLSG DTYEAVVTAVDPVADIATLR IQTKEPLPTLPLGRSADV RQ GEFV VAMGSPFALQNTITSG IVSSAQRPARDLGLPQTNVE YIQTDA AIDFGNSGGPLVNL DGEVIGVNTMKVTAGISFAI PSDRLREFLHRGEKKNSSSG IGSQRRYIGVMMLTLSPS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

