

## Human DJ-1/PARK7 protein (Recombinant) (No tag) (STJP000817)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mg/ml filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">11315</a>
Gene Symbol	<a href="#">PARK7</a>
UniProt ID	<a href="#">PARK7_HUMAN</a>
Immunogen Region	Met1-Asp189
Immunogen Sequence	MASKRALVILAKGAEEMETV IPVDVMRRAGIKVTVAGLAG KDPVQCSRDVVICPDASLED AKKEGPYDVVVLPGGNLGAQ NLSESAAVKEILKEQENRKG LIAAICAGPTALLAHEIGFG SKVTTHPLAKDKMMNGGHYT YSENRVEKDGLILTSRGPST SFEFALAIVEALNGKEVAAQ VKAPLVLKD

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

