

## Human MANF protein (Recombinant) (C-His) (STJP000819)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM PIPES, 500mM NaCl, pH 6.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7873</a>
Gene Symbol	<a href="#">MANF</a>
UniProt ID	<a href="#">MANF_HUMAN</a>
Immunogen Region	Leu25-Leu182
Immunogen Sequence	LRPGDCEVCISYLGRFYQDL KDRDVTFSPATIENELIKFC REARGKENRLCYYGATDDA ATKIINEVSKPLAHHIPVEK ICELKKKDSQICELKYDKQ IDLSTVDLKKLRVKELKKIL DDWGETCKGCAEKSDYIRKI NELMPKYAPKAASARTDL

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

