

Human FABP4/A-FABP/ALBP (A29T) protein (Recombinant) (No tag) (STJP000820)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2167
Gene Symbol	FABP4
UniProt ID	FABP4_HUMAN
Immunogen Region	Cys2-Ala132
Immunogen Sequence	CDAFVGTLVSSNFDDYM KEVGVGFTTRKVAGMAKPNM IISVNGDVITIKSESTFKNT EISFILGQEFDEVTADDRKV KSTITLDGGVLVHVQKWDGK STTIKRKREDDKLIVVECVMK GVTSTRVYERA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

