

Human Apolipoprotein A-I/APOA1 protein (Recombinant) (No tag) (STJP000822)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	335
Gene Symbol	APOA1
UniProt ID	APOA1_HUMAN
Immunogen Region	Asp25-Gln267
Immunogen Sequence	DEPPQSPWDRVKDLATVYVD VLKDSGRDYVSQFEGSALGK QLNKLLDNWDSVTSTFSKL REQLGPVTQEFWDNLEKETE GLRQEMSKDLEEVKAKVQPY LDDFQKKWQEEEMELYRQKVE PLRAELQEGARQKLHELQEK LSPLGEEMRDRARAHVDALR THLAPYSDELRRQLAARLEA LKENG GARLA EYHAKATEHL STLSEKAKPAEDLRQG LLP VLESFKVSFLSALEEYTKK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

