

Human NEDD8 protein (Recombinant) (C-His) (STJP000834)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4738
Gene Symbol	NEDD8
UniProt ID	NEDD8_HUMAN
Immunogen Region	Met1-Gly76
Immunogen Sequence	MLIKVKLTGKEIDIEPT DKVERIKERVEEKEGIPPQQ QRLIYSGKQMNDKTAADYK ILGGSVLHLVLALRGG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

