

Human ISG15 protein (Recombinant) (No tag) (STJP000836)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 50mM HEPES, 100mM NaCl, pH 7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9636
Gene Symbol	ISG15
UniProt ID	ISG15_HUMAN
Immunogen Region	Gly2-Gly157
Immunogen Sequence	GWDLTVKMLAGNEFQVSLSS SMSVSELKAQITQKIGVHAF QQRLAVHPSPGVALQDRVPLA SQGLPGSTVLLVVDKCDEP LSILVRNNKGRSSTYEVRLT QTV AHLKQQVSGLEGVQDDL FWLTFEGKPLEDQLPLGEYG LKPLSTVFMNLRRLGG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

