

Human FGF-12 protein (Recombinant) (C-His) (STJP000840)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, 1mM DTT, 300mM NaCl, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2257
Gene Symbol	FGF12
UniProt ID	FGF12_HUMAN
Immunogen Region	Met1-Thr181
Immunogen Sequence	MESKEPQLKGIVTRLFSQQG YFLQMHPDGTIDGTKDENS YTLFNLIPVGLRVVAIQGVK ASLYVAMNGEGYLYSSDVFT PECKFKESVFENYYVIYSST LYRQQESGRAWFLGLNKEGQ IMKGNRVKTKPSSHFPKP IEVCMYREPSLHEIGEKQGR SRKSSGTPTMNGGKVVNQDS T

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

