

Human R-spondin-1 protein (Active Recombinant) (C-His) (STJP000847)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	284654
Gene Symbol	RSPO1
UniProt ID	RSPO1_HUMAN
Immunogen Region	Arg31-Ala263
Immunogen Sequence	RISAEGSQACAKGCELCSEV NGCLKCSPKLFILLERNDIR QVGVCLPSCPPGYFDARNPD MNKCIKCKIEHCEACFSHNF CTKCKEGLYLHKGRCPYACP EGSSAANGTMECSSPAQCEM SEWSPWGPCSKKQQLCGFRR GSEERTRRVLHAPVGDHAAC SDTKETRRCTVRRVPCPEGQ KRRKGGQGRRENANRNLARK ESKEAGAGSRRRKGGQQQQQ QGTVGPLTSAGPA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

