

Human Fc-gamma RIIa/CD32a (H167R) protein (Recombinant) (C-His) (STJP000857)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2212
Gene Symbol	FCGR2A
UniProt ID	FCG2A_HUMAN
Immunogen Region	Ala36-Ile218
Immunogen Sequence	AAPPKAVLKLEPPWINVLQE DSVTLTCQGARSPESDSIQW FHNGNLIPTHTQPSYRFKAN NNDSGEYTCQTGQTSLSDPV HLTVLSEWLVLQTPHLEFQE GETIMLRCHSWKDKPLVKVT FFQNGKSQKFSHLDPTFSIP QANHSMSGDYHCTGNIGYTL FSSKPVTITVQVPSMGSSSP MGI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

