

## Human VEGFR-2/KDR/CD309 protein (Recombinant) (C-His) (STJP000860)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3791</a>
Gene Symbol	<a href="#">KDR</a>
UniProt ID	<a href="#">VGFR2_HUMAN</a>
Immunogen Region	Ala20-Glu764
Immunogen Sequence	ASVGLPSVSLDLPRLSIQKD ILTIKANTTLQITCRGQRDL DWLWPNQSGSEQRVEVTEC SDGLFCKTLTIPKVIQNDTG AYKCFYRETDLASVIYVYVQ DYRSPFIASVSDQHGVVYIT ENKNKTVVIPCLGSISNLNV SLCARYPEKRFVDPGNRISW DSKKGFTIPSYMISYAGMVF CEAKINDESYSQSIMYIVVVV GYRIYDVVLSPSHGIELSVG EKLVLNCTARTELVGIDF

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

