

Human Tyrosine-protein kinase receptor UFO/Axl protein (Recombinant) (C-His) (STJP000863)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	558
Gene Symbol	AXL
UniProt ID	UFO_HUMAN
Immunogen Region	Glu33-Pro449
Immunogen Sequence	EESPFVGNPGNITGARGLTG TLRQCLVQGEPEVHWLRD GQILELADSTQTQVPLGEDE QDDWIVVSQLRITSLQLSDT GQYQCLVFLGHQTFVSPQGY VGLEGLPYFLEEPEDRTVAA NTPFNLSCQAQGPPEPVDLL WLQDAVPLATAPGHGPQRSL HVPGLNKTSSFSCAEHNAKG VTTSRATITVLPQQPRNLH LVSQPTELEVAWTPGLSGI YPLTHCTLQAVLSDDGMGI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

