

Human MBP-C/MBL-2 protein (Recombinant) (C-His) (STJP000871)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4153
Gene Symbol	MBL2
UniProt ID	MBL2_HUMAN
Immunogen Region	Glu21-Ile248
Immunogen Sequence	ETVTCEDAQKTCPAVIACSS PGINGFPGKDGRDGTKEKG EPGQGLRGLQGPPGKLGPPG NPGPSGSPGPKGQKGDGPKS PDGDSSLAASERKALQTEMA RIKKWLTFSLGKQVGNKFFL TNGEIMTFEKVKALCVKFQA SVATPRNAAENGAIQNLIKE EAFLGITDEKTEGQFVDLTG NRLTYTNWNEGEPNAGSDE DCVLLKNGQWNVPCSTSH LAVCEFFI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

