

Human IL-18R1/CD218a protein (Recombinant) (C-hFc&His) (STJP000873)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8809
Gene Symbol	IL18R1
UniProt ID	IL18R_HUMAN
Immunogen Region	Glu20-Arg329
Immunogen Sequence	ESCTSRPHITVVEGEPFYLK HCSCSLAHEIETTTKSWYKS SGSQEHVELNPRSSRIALH DCVLEFWPVELNDTGSYFFQ MKNYTKQWKLNVIRRNKHSC FTERQVTSKIVEVKKFFQIT CENSYYQTLVNSTSLYKNCK KLLLENNKNPTIKKNAEFED QGYYSVCVHFLHHNGKLFNIT KTFNITIVEDRSNIVPVLGG PKLNHVAVELGKNVRLNCSA LLNEEDVIYWMFGEENGSD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

