

Human MCAM/CD146 protein (Recombinant) (C-His) (STJP000877)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4162
Gene Symbol	MCAM
UniProt ID	MUC18_HUMAN
Immunogen Region	Val24-Gly559
Immunogen Sequence	VPGEAEQPAPLVEVEVGST ALLKCGLSQSQGNLSHVDWF SVHKEKRTLIFRVRQQGQGS EPGEYEQRLSLQDRGATLAL TQVTPQDERIFLCQGKRPRS QEYRIQLRVYKAPEEPNIQV NPLGIPVNSKEPEEVATCVG RNYPIPVQVIWYKNGRPLKE EKNRVHIQSSQTVESSGLYT LQSILKAQLVKEDKDAQFYC ELNYRLPSGNHMKESREVTV PVFYPTKVVLEVEPVGML

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

