

Human MMP-9 protein (Recombinant) (C-His) (STJP000879)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4318
Gene Symbol	MMP9
UniProt ID	MMP9_HUMAN
Immunogen Region	Ala20-Asp707 Q279R
Immunogen Sequence	APRQRQSTLVLFPGDLRTNL TDRQLAEEYLYRYGYTRVAE MRGESKSLGPALLLLQKQLS LPETGELDSATLKAMRTPRC GVPDLGRFQTFEGDLKWHHH NITYWQNYSEDLPRVIDD AFARAFALWSAVTPLTFTRV YSRDADIVIQFGVAEHGDGY PFDGKDGLLAHAFPPGPGIQ GDAHFDDDELWSLGKGVVVP TRFGNADGAACHFPFIFEGR SYSACTTDGRSDGLPWCST

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

