

Human Basigin/CD147 protein (Recombinant) (C-hFc&His) (STJP000880)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	682
Gene Symbol	BSG
UniProt ID	BASI_HUMAN
Immunogen Region	Thr25-His205
Immunogen Sequence	TVFTTVEDLGSKILLTCSLN DSATEVTGHRWLKGGVVLKE DALPGQKTEFKVDSDDQWGE YSCVFLPEPMGTANIQLHGP PRVKAVKSSEHINEGETAML VCKSESVPPVTDWAWYKITD SEDKALMNGSERFFVSSSQ GRSELHIENLNMEADPGQYR CNGTSSKSGSDQAITLRVRS H

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

