

Human SLAMF2/CD48 protein (Recombinant) (C-hFc&His) (STJP000883)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	962
Gene Symbol	CD48
UniProt ID	CD48_HUMAN
Immunogen Region	Gln27-Ser220
Immunogen Sequence	QGHVHMTVVSGSNVTLNIS ESLPENYKQLTWFYTFDQKI VEWDSRKSKYFESKFKGRVR LDPQSGALYISKVQKEDNST YIMRVLKKTGNEQEWKIKLQ VLDVPVKPVIKIEKIEDMDD NCYLKLSCVIPGESVNYTWY GDKRPFPKELQNSVLETTLM PHNYSRCYTCQVNSNSVSSKN GTVCLSPCTLARS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

