

Human CDH17/LI-cadherin protein (Recombinant) (C-His) (STJP000884)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1015
Gene Symbol	CDH17
UniProt ID	CAD17_HUMAN
Immunogen Region	Gln23-Met787
Immunogen Sequence	QEGKFSGPLKPMTFSIYEGQ EPSQIIFQFKANPPAVTFEL TGETDNIFVIEREGLLYNR ALDRETRSTHNLQVAALDAN GIIVEGVPVITIKVKDINDN RPTFLQSKYEGSVRQNSRPG KPFLYVNATDLDDPATPNGQ LYYQIVQLPMINNVMYFQI NNKTGAISLTREGSQELNPA KNPSYNLVISVKDMGGQSEN SFSDTTSVDIIVTENIWKAP KPVEMVENSTDPHIKITQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

