

Human Clusterin/Apo-J/CLU protein (Recombinant) (C-His) (STJP000885)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1191
Gene Symbol	CLU
UniProt ID	CLUS_HUMAN
Immunogen Region	Asp23-Arg227
Immunogen Sequence	DQTVSDNELQEMSNQGSKYV NKEIQNAVNGVKQIKTLIEK TNEERKTLLSNLEEAKKKKE DALNETRESETKLKELPGVC NETMMALWEECKPCLKQTCM KFYARVCRSGSGLVGRQLEE FLNQSSPFYFWMNGDRIDSL LENDRQQTHMLDVMQDHFSR ASSIIDELFQDRFFFTREPQD TYHYLPFSLPHRRPHFFFPK SRIVRSLMPFSPYEPLNFHA MFQPFLEMIHEAQQAMDIH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

