

Human Coagulation factor III/CD142 protein (Recombinant) (C-His) (STJP000888)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2152
Gene Symbol	F3
UniProt ID	TF_HUMAN
Immunogen Region	Gly34-Glu251
Immunogen Sequence	GTTNTVAAYNLTWKSTNFKT ILEWEPKPVNQVYTVQISTK SGDWKSCKFYTTDTECDLTD EIVKDVKQTYLARVFSYPAG NVESTGSAGEPLYENSPEFT PYLETNLGGPTIQSFEQVGT KVNVTVEDERTLVRNNNTFL SLRDVFGKDLIYTYWYKSS SSGKKTAKTNTNEFLIDVDK GENYCFVQAVIPSRVTNRK STDSPVECMGQEKGEFRE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

