

Human DDR2/CD167b protein (Recombinant) (C-His) (STJP000891)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4921
Gene Symbol	DDR2
UniProt ID	DDR2_HUMAN
Immunogen Region	Gln24-Arg399
Immunogen Sequence	QVNPACRYPLGMSGGQIPD EDITASSQWSESTAAKYGRLL DSEEGDGAWCPEIPVEPDDL KEFLQIDLHLHFITLVGTQ GRHAGGHGIEFAPMYKINYS RDGTRWISWRNRHGKQVLDG NSNPYDIFLKDLEPPIVARF VRFIPVTDHSMNVCMRVELY GCVWLDGLVSYNAPAGQQFV LPGGSIYLNDSVYDGAVGY SMTEGLGQLTDGVSGLDDFT QTHEYHVVWPGYDYYGWRNE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

