

## Human IL-3RA/CD123 protein (Recombinant) (C-His) (STJP000897)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3563</a>
Gene Symbol	<a href="#">IL3RA</a>
UniProt ID	<a href="#">IL3RA_HUMAN</a>
Immunogen Region	Lys20-Arg305
Immunogen Sequence	KEDPNPPITNLRMKAKAQQQL TWDLNRNVTDIECVKDADYS MPAVNNSYCQFGAISLCEVT NYTVRVANPPFSTWILFPEN SGKWPWAGAENLTCWIHDVDF LSCSWAVGPGAPADVQYDLY LNVANRRQYECLEHYKTDQA GTRIGCRFDDISRLSSGSQS SHILVRGRSAAFIPCTDKF VVFSQLIETPPNMTAKCNK THSFMHWKMRSHFNRFKFRYE LQIQKRMQPVITEQVRDRT

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

