

Human CEACAM3/CD66d protein (Recombinant) (C-His) (STJP000907)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1084
Gene Symbol	CEACAM3
UniProt ID	CEAM3_HUMAN
Immunogen Region	Met1-Gly155
Immunogen Sequence	MGPPSASPHRECIPWQGLLL TASLLNFWNPTTAKLTIES MPLSVAEGKEVLLLHNLVLPQ HLFGYSWYKGERVDGNLIV GYVIGTQQATPGAAYSGRET IYTNASLLIQNVLTQNDIGFY TLQVIKSDLVNEEATGQFHV YQENAPGLPVGAVAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

