

Human FGFR-5/FGFRL1 protein (Recombinant) (C-His) (STJP000909)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53834
Gene Symbol	FGFRL1
UniProt ID	FGRL1_HUMAN
Immunogen Region	Ala 25-Pro 378
Immunogen Sequence	ARGPPKMADKVVPRQVARLG RTVRLQCPVEGDPPLTMWT KDGRTHSGWSRFRVLPQGL KVKQVEREDAGVYVCKATNG FGSLSVNYTLVVLDDISPGK ESLGPDSSSGGQEDPASQQW ARPRFTQPSKMRRRVIARPV GSSVRLKCVASGHPRPDITW MKDDQALTRPEAAEPRKKKW TSLSLKNLRPEDSGKYTCRVS NRAGAINATYKVDVIQRTRS KPVLTGTHPVNTTVDFGGT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

