

Human LIMP II/SCARB2/CD36L2 protein (Recombinant) (C-His) (STJP000910)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	950
Gene Symbol	SCARB2
UniProt ID	SCRB2_HUMAN
Immunogen Region	Arg27-Thr432
Immunogen Sequence	RVFQKAVDQSIIEKIVLRNG TEAFDSWEKPLPVYTQFYF FNVTNPEEILRGETPRVEEV GPYTYRELRNKANIQFGDNG TTISAVSNKAYVFERDQSVG DPKIDLIRTLNIPVLTVIEW SQVHFLREIIEAMLKAYQQK LFTVHTVDELLWGYKDEILS LIHVFRPDISPYFGLFYEKN GTNDGDYVFLTGEDSYLNFT KIVEWNGKTSLDWWITDKCN MINGTDGDSFHPLITKDEV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

