

Human Glypican-3/GPC3 protein (Recombinant) (C-His) (STJP000913)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2719
Gene Symbol	GPC3
UniProt ID	GPC3_HUMAN
Immunogen Region	Gln25-His559
Immunogen Sequence	QPPPPPPDATCHQVRSFFQR LQPLKWWPETPVPGSDLQV CLPKGPTCCSRKMEEKYQLT ARLNMEQLLQSASMEKFLI IQNAAVFQEAFFIIVRHAKN YTNAMFKNNYPSLTPQAFEF VGEFFTDVSLYILGSDINVD DMVNELFDSLFPVIYTQLMN PGLPDSALDINECLRGARRD LKVFNGFPKLIMTQVSKSLQ VTRIFLQALNLGIEVINTTD HLKFSKDCGRMLTRMWYCS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

