

Human Complement C1r subcomponent protein (Recombinant) (C-His) (STJP000917)

GENERAL INFORMATION

| | |
|---------------|--------------|
| Product Type | Proteins |
| Host / Source | HEK293 cells |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 715 |
| Gene Symbol | C1R |
| UniProt ID | C1R_HUMAN |
| Immunogen Region | Met1-Asp705 |
| Immunogen Sequence | MWLLYLLVPALFCRAGGSIP IPQKLFGEVTSPLFPKPYPN NFETTTVITVPTGYRVKLVF QQFDLEPSEGCFYDYVKISA DKKSLGRFCGQLGSPLGNPP GKKEFMSQGNKMLLTFHTDF SNEENGTIMFYKGFLLAYYQA VDLDECASRSKLGEEDPQPQ CQHLCHNYVGGYFCSCRPGY ELQEDRHSCQAECSSSELYTE ASGYISSLEYPRSYPPDLRC NYSIRVERGLTLHLKFLEP |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|

