

Human FGFR-3/CD333 protein (Recombinant) (C-hFc&His) (STJP000919)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2261
Gene Symbol	FGFR3
UniProt ID	FGFR3_HUMAN
Immunogen Region	Glu23-Gly375
Immunogen Sequence	ESLGTEQRVVGRAAEVPGPE PGQQEQLVFGSGDAVELSCLP PPGGGPMGPTVWVKDGTGLV PSEKLVGPGQRLQVLNASHE DSGAYSCRQLTQRVLCHFS VRVTDAPSSGDDDEDGEDEAE DTGVDTGAPYWTRPERMDKK LLAVPAANTVRFRCPAAGNP TPSISWLKNGREFRGEHRIG GIKLRHQWLSLVMESVPSD RGNVTCVVENKFGSIRQTYT LDVLEKSPHRPILQAGLPA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

